Your Endo Specialist™

HEROfill®

Your Endo Specialist™

Operating instructions

HEROfill® Obturator and Verifier Operating Instructions

Classification

The HEROfill® Obturator is a Class IIa medical device, in compliance with MDD (Medical Devices Directive 93/42/EEC). The HEROfill® Obturator must only be used by dentists.

Read and follow the instructions carefully!

General precautions

- The product should only be used by dentists.
- The product should only be used for root canal procedures in teeth.
- A deformed or defective product should never be used, but should be discarded.
- The manufacturer disclaims all responsibility and liability for injury or damage to persons or property caused by faulty and/or inappropriate use. For information on the correct use of the product, please see the step-by-step guide (section 2).
- The product should be stored in the original packaging.

Warnings - HEROfill® Obturator

- After heating, the Gutta Percha is hot. Avoid any direct contact between the HEROfill® Obturator and skin or mucous membranes. The HEROfill® Obturator should be inserted directly into the root canal.
- It is always recommended to take a control x-ray. If the x-ray is not satisfactory, take another x-ray from a different angle. If the filling still cannot be satisfactorily documented, remove the HEROfill® Obturator.

1. Description

A HEROfill® Obturator is a carrier of biocompatible plastic, which is coated with thermoplastic Gutta Percha.

The sizes are 20 to 45 according to ISO standard. The HEROfill® Obturator sizes correspond to the ISO standards of root canal files. A single HEROfill® Obturator is all that is needed to obturate a root canal fully.

The result should always be a perfectly placed root canal filling with a tight apical seal.

A HEROfill® Obturator consists of three parts:

1. Plastic handle with a metal insertion pin.

The plastic handle is colour-coded according to the size of the HEROfill® Obturator. The stainless steel insertion pin is 9 mm long, and is permanently attached to the plastic handle. The metal insertion pin is placed 6 mm down in the plastic carrier. When the metal insertion pin is removed it will leave a hollow part of 6 mm in the plastic carrier.

2. Outer layer of thermoplastic Natural Gutta Percha.

The Gutta Percha is reversible thermoplastic. This means that it becomes soft and highly adhesive when heated to temperatures above approx. 100°C. When cooled, it returns to a firm, rigid state. Heating up will once again make the Gutta Percha soft and adhesive.

3. Plastic carrier with endo stop.

The tapered plastic carrier is sufficiently flexible to negotiate curved canals.

Your Endo Specialist™

The HEROfill® Obturator plastic carrier is slightly tapered. This ensures adequate back flow of the heated Gutta Percha. The length of the carrier is 24 mm. The coronal part has length indication for each mm (18-23 mm).

The coronal 6 mm of the carrier is hollow to accommodate the metal insertion pin. A rubber endo stop is provided to indicate the working distance measurement. The plastic carrier and Gutta Percha are radiopaque.

2. A step-by-step guide

After having ideally used the One Shape® or Revo-S™ methods for shaping and before using the HEROfill® Obturator, assure that you have achieved the following:

- A properly cleaned and shaped canal – the key to any successful root canal filling.
- A dry root canal.
- An established working length. The plastic carrier is flexible enough to negotiate curved canals. It is highly recommended that rotary files with a higher conicity are utilized to create a sufficient taper for the HEROfill® Obturator. When using files with a conicity of 2%, it is advised to combine them with a procedure using Gates Glidden Drills, orifice openers (ENDOFLARE® for example) or similar files which can be used to create a sufficient opening in the coronal part of the canal. The use of HEROfill® Verifiers is recommended.

STEP 1: Determine the working length and choose the proper HEROfill® Obturator size

The working length is measured all the way to the apex (A1). We recommend measuring the working length from the edge of the cusp. The working length may be determined by use of x-ray, file and electronic apex locator. Generally, the appropriate size HEROfill® Obturator will be the same size as the last file used at the apex of the canal. *Tip: For very narrow and/or highly calcified canals, it may be useful to select one size smaller than the last file used at apex.*

If you do not use rotary files with a taper of 4% or more, you should always verify the size of the cleaned root canal with a HEROfill® Verifier before inserting the HEROfill® Obturator.

Mark your working length on your HEROfill® Obturator. When the HEROfill® Verifier has a slightly loose fit in the apical third, the rubber stop is placed at the working length. Afterwards, the working length is transferred to the HEROfill® Obturator.

STEP 2: Heating of the HEROfill® Obturator

Place the selected HEROfill® Obturator in the HEROfill® Oven and activate it (A2, A3). For more information on how to use the HEROfill® Oven, please consult the HEROfill® Oven Operating Instructions.

STEP 3: Application of sealer

While the HEROfill® Obturator is heating, mix any heat resistant sealer (usually a non-eugenol sealer – MM-Seal™ recommended) with long working time. Make sure that the root canal is completely dry before applying the sealer.

Place a very thin coat of sealer on the wall of the canal. Too much sealer is neither desirable nor necessary. Utilize a Paper Point (MM Paper Points recommended) or a HEROfill® Verifier to apply a thin layer of sealer on the canal walls. When obturating multiple canals, apply sealer in all canals at once.

STEP 4: Insertion of the HEROfill® Obturator into the canal

When the heater indicates that the HEROfill® Obturator is ready for use, the HEROfill® Obturator is removed from the HEROfill® Oven. For further information on how the HEROfill® Oven works, please consult the appropriate heater manual. Without twisting the handle, immediately insert the HEROfill® Obturator into the canal to the working distance using firm and steady pressure (A4).

Your Endo Specialist™

STEP 5: Confirmation X-ray

Confirm your root canal treatment by taking X-rays.

STEP 6: Removal of excess central carrier and Gutta Percha

When the Gutta Percha has hardened (after 2-3 min.) remove the handle in one of two ways (A5).

1. Lean the handle and metal insertion pin sharply to one side to break off the excess. You may pull the handle up as needed and break for each millimetre.
2. Stabilize the plastic carrier with a finger, and cut away the handle and excess plastic core with a small inverted cone bur.

Be careful to prevent vertical movement of the carrier to minimize the chance of disturbing the apical seal. Trim away the extra Gutta Percha using an endo excavator or similar instrument.

3. Technique tips

It is very easy to use the HEROfill® system. Feedback from experienced HEROfill® system users, however, will hopefully make it even easier for you to get started.

Obturing multi-rooted cases

There will almost always be more than enough Gutta Percha on each HEROfill® Obturator. When obturating the canal, the excess Gutta Percha and sealer will backfill to the coronal part of the tooth, possibly blocking the adjacent opening. However, you can avoid this by:

- Obturating the shortest canal first.
- Placing Paper Points or HEROfill® Verifier in the other canals prior to obturation to prevent the following root canals being blocked by Gutta Percha. Remove Paper Points or size HEROfill® Verifiers as each canal is obturated.
- While there is sufficient Gutta Percha on each HEROfill® Obturator to fill even extremely long canals, wide canals or internally resorbed canals, it may be too much in some cases. If it is obvious that there will be an excess of Gutta Percha to fill the canal, use a sharp blade to trim some of the material from the coronal end of the central carrier prior to heating.
- In some cases it can be beneficial to remove the handle/metal insertion pin from the actual plastic carrier, especially in multi-rooted teeth where the extra space will make it easier to place the following HEROfill® Obturators in the second and third canal.

Prior to heating of the HEROfill® Obturator(s) and prior to marking the working length, always assure that the handle/metal insertion pin can be separated from the plastic carrier. This is done by twisting the handle/metal insertion pin out of the plastic carrier. Re-insert the handle/metal insertion pin into the plastic carrier.

Mark your working length and proceed to STEP 2 – STEP 4 (section 2).

Stabilize the plastic carrier with a cotton plier while the handle/metal insertion pin is twisted out. Any excess material is removed as described in STEP 5 (section 2).

Use of sealer

As with any other root filling technique, when using the HEROfill® System, always apply sealer in the canal first. However, there is a significant difference when using HEROfill® System. Only apply a very small amount using a HEROfill® Verifier or a Paper Point. A root spiral should not be utilized when using HEROfill® Obturators as there is a tendency to add too much sealer to the canal. Remember; only apply a very thin layer of sealer on the canal walls.

4. Post space preparation

Post space creation can be accomplished by removing the coronal portion of the obturation. Since 6 mm of the plastic core is hollow, this procedure should be easier compared to solid core obturators. Smooth round burs, preferably a ball burnisher with no cutting flutes, are the most effective since they do not cut, but instead “melt” the gutta

Your Endo Specialist™

percha and central core by friction. Cutting burs may be used, but care should be taken to prevent perforation.

Using a round bur, or ball burnisher:

1. Use a high speed hand piece with water coolant.
2. Begin with the largest size round bur that will fit in the access opening of the canal without touching the walls.
3. Place the bur in the centre of the hollow core using gentle, relaxed pressure. Intermittently pulse the foot controller to allow the bur to remove or melt the plastic core.
4. Change to smaller diameter burs as the canal tapers to prevent damage to the walls until the desired depth is reached.

5. Retreatment of an obturated root canal

Removal of the HEROfill® Obturator is relatively simple prior to the Gutta Percha cooling.

Simply grasp the excess carrier with cotton pliers and remove it.

Once the Gutta Percha has cooled, or for retreatment at a later appointment, removal is easiest achieved by using a Peeso bur, Gates or similar:

1. Place the bur between the canal wall and the plastic carrier, so that the rotational direction will be counter-clockwise. Use the instrument at very low speeds (250-400 rpm), and work slowly around the plastic carrier until you feel resistance. The instrument will grasp between the plastic carrier and the canal wall. The plastic carrier will be coiled out of the canal after a short time.
2. If necessary, a small Hedstrom hand file can be worked down apically along the plastic carrier.
3. When the plastic carrier has been pulled out of the canal, the Gutta Percha is removed by traditional means.

6. Packaging

HEROfill® products are available as follows:

- The HEROfill® Obturator is packed individually in a sealed plastic box. One box consists of 1 x 6 pcs.
- The HEROfill® Verifiers are packed in the same way as the HEROfill® Obturator.

7. Labelling

	Manufacturer
	See instruction manual
	Single use
	Expiry date
	Lot number
	CE Class I HEROfill® Verifiers
	CE 2460 Class IIa HEROfill® Obturators

Contact your local supplier for more information.

Your Endo Specialist™

HEROfill® Oven Operating Instructions

Classification

The HEROfill® Oven is a Class I medical device, in compliance with MDD (Medical Devices Directive 93/42/EEC) Section IX.

The heater must only be used in dental clinics by dental personnel, and only for heating the HEROfill® Obturator. The product is classified according to DS EN 61010-1-2010 and therefore must not be used within reach of patients.

The HEROfill® Oven has a heat chamber that consists of 5 printed circuit boards with copper lanes. By applying power to the copper lanes heat is generated.

The oven is driven by an external power supply.

The HEROfill® Oven must only be serviced and repaired by MICRO-MEGA® or an approved dealer otherwise the warranty will be void.

Warning

Do not touch the inside of the heat chamber before the oven is completely cooled.

The Heater must be used for its intended purpose and is protected from overheating by two independent circuits.

Technical specifications

POWER SUPPLY	
Input power	100-240V AC, 50/60Hz
Output power	9V DC, 1.7A
SAFETY	
Auto turn off	The heater turns off automatically after 9 hours

(Normal) Environmental conditions

This standard applies to equipment designed to be safe at least under the following conditions:

a) indoor use;
b) altitude up to 2000 m;
c) temperature 5°C to 40°C;
d) maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% humidity at 40°C;
e) Mains supply fluctuations + 10% of normal voltage;
f) Overvoltage Category II;
g) Temporary overvoltages occurring on the mains supply;
h) Pollution degree 2

Heating times

HEROfill® Obturator heating time is 60 seconds.

Preheating of the heater takes 3-4 minutes depending on surrounding temperature.

1. Start using the heater

1. Connect the power supply to the DC inlet of the heater and the mains and switch on power (B1).
2. Indicator A emits a red light while the HEROfill® Oven preheats (B2). When the HEROfill® Oven is ready, the indicator emits a green light, and you will hear a beep.

Notice: The HEROfill® Oven can be left on the entire day.

Your Endo Specialist™

2. Use

When the heater is ready for use, up to four HEROfill® Obturators at a time may be placed in the appropriate slots. Place the HEROfill® Obturators in slots marked 1 to 4 (B3). *Tip: Place the HEROfill® Obturator to be used first in slot #1, the second in #2 etc. - generally you should begin with the smallest size HEROfill® Obturator.*

1. Activate the timer by pressing the **ON** icon (B4).
2. Indicator **B** lights up red during the heating cycle.
3. When ready, the indicator changes to green and you will hear a beep (B4).

Let the HEROfill® Obturator stay in the heater until you are ready to insert it (max. 15 min).

4. It is possible to turn off the HEROfill® Oven by holding down the **ON** button for 2 seconds (B4). The HEROfill® Oven beeps three times to indicate it is switched off. Hold down the **ON** button for 2 seconds to reactivate the heater.

Tip: You do not have to turn off the HEROfill® Oven after finished treatment; the HEROfill® Oven can be left on standby all day.

Indicators - light and sound

	Time	Heat	Sound
Preheating	●	●	-
Ready	●	●	Long beep
Activation	●	●	Short beep
Obturator ready	●	●	Long beep
Turn off	●	●	3 x short beep
Power on	●	● / ●	Long beep
Auto turn off	-	-	Long beep

3. Cleaning

When the heater has cooled completely, gutta percha residue may easily be removed. The top of the heater may be removed for cleaning by rotating the top counter-clockwise 1 until the top is loose and can be removed 2 (B5).

The top may be washed in a dishwasher and/or autoclaved. Gutta percha residue may be removed from the heat chamber by the use of a blunt plastic instrument.

The main part of the heater may be wiped with a slightly damp cloth, using standard cleaning solutions and surface disinfection products.

Sterilization guide for removable top part


WARNINGS	DISCHARGE DAMAGED PRODUCTS
Limitations on reprocessing	250 cycles
Instructions	
Point of use	Remove excess soil with disposable cloth/paper tissue.
Preparation for decontamination	No specific requirements
Cleaning: Automated	Not recommended
Cleaning: Manual	Equipment: Detergent, brush, distilled or demineralised water.

Your Endo Specialist™

	<ol style="list-style-type: none"> 1. Rinse excess soil from top 2. Using brush, apply detergent solution to all surfaces 3. Rinse under clear running water
Desinfection	Not recommended
Drying	Do not exceed 134° C
Maintenance, inspection and testing	Discharge damaged products
Packaging	Only pack individually
Sterilization	Vacuum autoclave, minimum 18 minutes at 134° C, pressure 2.2 bar. Do not exceed 134° C.
Storage	No particular requirements

4. Labelling

	Manufacturer
	See instruction manual
	Serial number
	Medical Electrical Equipment Directive, DS/EN 60601-1:2006 Electric shock protection, Type B
	Medical Devices Directive, MDD Directive 93/42/EEC Medical device, Class I
	Canadian Standards Association, CSA Approved for use in Canada and the USA
	Waste Electrical and Electronic Equipment Directive, WEEE Directive 2002/96/EC Please return the product for recycling

Accessories

HEROfill® Obturator References

- HEROfill® Obturator n°20 - Ref. 20138501 N
- HEROfill® Obturator n°25 - Ref. 20138502 N
- HEROfill® Obturator n°30 - Ref. 20138503 N
- HEROfill® Obturator n°35 - Ref. 20138504 N
- HEROfill® Obturator n°40 - Ref. 20138505 N
- HEROfill® Obturator n°45 - Ref. 20138506 N

Accessories References

- HEROfill® Verifier n°20 - Ref. 20139501 N
- HEROfill® Verifier n°25 - Ref. 20139502 N
- HEROfill® Verifier n°30 - Ref. 20139503 N
- HEROfill® Verifier n°35 - Ref. 20139504 N
- HEROfill® Verifier n°40 - Ref. 20139505 N
- HEROfill® Verifier n°45 - Ref. 20139506 N

Your Endo Specialist™

Service and repair

Inspections and repairs may only be carried out by an approved dealer or by MICRO-MEGA®. Your MICRO-MEGA® product is an accurate medical device, produced under maximum quality and test conditions.

To guarantee the personal follow-up and traceability of medical devices, please register the guarantee of your HEROfill® on our website (www.micro-mega.com) within 10 days of your purchase. This return ensures that the after-sales procedures for your device will be activated so that MICRO-MEGA® can give you full satisfaction with the service you expect and manage your warranty correctly.

MICRO-MEGA® guarantees your HEROfill® Oven for two years against any hidden manufacturing defect, excluding faults due to:

- Improper use.
- Inappropriate maintenance or application, not in accordance with our recommendations and instructions for use.
- Accidental damage (fall, impact, etc.).
- Attempted intervention or repair by personnel not approved by MICRO-MEGA®.
- Attempted product modification.

Waste handling within the EC

In compliance with the WEEE Directive.

Waste disposal within the EC:

1. WEEE (Waste of Electrical and Electronic Equipment) cannot be disposed with normal household waste.
2. WEEE must be collected separately.
3. The equipment has to be treated at an approved re-treatment centre.
4. The equipment cannot be reused after its lifetime has ended.
5. WEEE pollutes the environment and thereby damages the health of living beings.

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