# STAT/S

Pièce à main haute vitesse / Turbines  
High-Speed Handpieces  
Turbinen  
Manipoli ad alta velocità  
Piezas de mano de alta velocidad

**Manuel de l'opérateur**  
**Operating instructions**  
**Benutzerhandbuch**  
**Manuale per operatori**  
**Manual del Operador**

**STAT/S SL 101.1**  
**ML 201.1**

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## Table of Contents

1	Important Information .....	18
	Applied Symbols.....	18
	Intended Use - Purpose .....	20
2	Accessories .....	22
3	Reprocessing Methods.....	22
4	Installation and Commissioning.....	26
5	Operation .....	27
6	Maintenance and Servicing .....	30
7	Specifications.....	31
8	Guarantee .....	32

## 1 Important Information

Dear user,  
We hope that you will enjoy using your new quality product. In order to ensure that you can work in a trouble-free, economical and safe manner, please comply with the instructions below.

### Applied Symbols

Appearing in the manual:

	Situations where failure to follow the instructions may lead to danger, damage to material or operating faults.
	Important information for operator and engineer.
	Sterilizable up to 135°C (275°F)
	CE mark (Communauté Européenne)

Appearing on the packaging:

	Fragile
	Keep dry
	Stacking restrictions
	Temperature range
	Air pressure
	Humidity
	Quantity

Modifications and improvements to the product as a result of new technical developments are possible.  
This does not imply any right to retrofitting of existing units.

## 1 Important Information

The manufacturer assumes no responsibility for damage arising through:

- External influences (poor quality of the media or faulty installation)
- Use of incorrect information
- Improper use
- Improperly performed repairs.

Repair and maintenance work - apart from the activities described in this operator's manual - may be performed only by qualified technical staff.

- In the event of modifications by third parties, the licences become null and void.
- Only use original parts and spare parts.

### **Endangerment to the practitioner and patient**

You must stop using the device immediately in the event of damage, irregular running noises, excessive vibrations, unusual generation of heat or if the drill bit is not firmly gripped.

- The turbine must be cleaned, serviced and stored dry if it is to be put into extended storage.

### **Risk due to incorrectly stored turbine**

Damage to the turbine from dropping

- After treatment, place the turbine in the hand-piece holder of the delivery unit.

Any claim under warranty shall be excluded if defects or the consequences thereof are due to manipulation or modifications to the product by the customer or by any third parties not authorized by MICRO-MEGA®.

## 1 Important Information

### Intended Use - Purpose

This turbine is:

- Only intended for dental treatment in the practice of dentistry. Any other type of use or alteration to the product is impermissible and can be hazardous.
- A turbine according to relevant national statutory regulations
- A dental contra-angle handpiece, air-driven, with light, push-button chuck, 1 - 4 internal water sprays
- Not approved for use in areas with an increased risk of explosion

According to these provisions, the turbine is only to be used by an experienced user for the described application in accordance with:

- The applicable health and safety regulations
- The applicable accident prevention regulations
- This operator's manual

According to these regulations, the user is required to:

- Use only tools which are free of faults and works properly
- Use only tools for the proper purpose
- Protect himself, the patient and third parties from danger
- Avoid contamination from the product

This turbine is intended for following applications:

- **Indications:** Intended for removal of carious material, reducing hard tooth structure, cavity preparation, finishing tooth preparations, restorations and for polishing teeth.
- **Contra indications:** not known
- **Adverse effects:** not known

## 1 Important Information

### General warnings

- To be used by dental practitioners only.
- Perform maintenance on a regular basis.
- Use good quality tools.
- It is recommended that a dental dam be used during treatment.
- Never dismantle the turbine. The device may only be dismantled by the approved distributor or at the factory.
- Never press the push button while the device is rotating.
- It is recommended that suitable protective gloves and eye protection be worn when handling the turbine and instruments.
- For extended periods of use, relubricate every 10 minutes.
- Contact the MICRO-MEGA® after sales department straight away in the event of a problem with the turbine.
- Always remove the diamond tip if any part of the diamond is missing.
- Check the instrument is present before starting up the turbine.

### Warnings regarding reprocessing (cleaning - disinfection - lubrication - sterilisation).

- Ensure that the turbine is sterilised and lubricated prior to initial use and between patients.
  - Remove the rotating parts (burs, files, etc.) before reprocessing.
  - Avoid soaking the turbine in a tub or ultrasonic cleaning system.
  - Avoid soaking the turbine in a special purpose detergent for longer than 15 minutes.
  - Avoid working with a turbine if it is overheating.
  - Use a detergent-disinfectant designed specifically for pre-disinfection of dynamic dental instruments, which contains neither sodium oxide nor aldehyde and which is aluminium compatible. Do not use detergents or solutions which are alkaline (pH > 7).
  - Use the concentration recommended by the pre-disinfection product manufacturer.
  - Observe the immersion time recommended by the pre-disinfection product manufacturer.
- Remove the turbine from the sterilisation bag immediately before use.

## 1 Important Information

Do not use vinegar to clean the spray nozzles.

If the turbine is to remain unused for an extended period of time, it must be cleaned and lubricated in accordance with the instructions and stored in a dry place. Stop using the turbine if it becomes damaged or overheats, if it makes an abnormal sound during use, if there is excessive vibration, if the drill becomes excessively out-of-round or does not fit properly. To avoid any further damage, contact the MICRO-MEGA® customer services department or your approved MICRO-MEGA® stockist.

## 2 Accessoires

Lubricating:	Ref :
Oil spray «500l»	50700008
Spray nozzle	50700032

## 3 Reprocessing Methods

### Reprocessing limitations

- The sterilisation process does not affect the performance of dynamic instruments. The product's life span normally depends on damage incurred during use.

### Place of use

- After use, remove the excess debris with a disinfectant wipe.

### Containment and transport

- It is advisable to reprocess the turbine as soon as possible after use or to immerse it in water to prevent debris from drying and to prevent blood and saliva from coagulating.

### 3 Reprocessing Methods

#### Manual cleaning

##### *(Steps to take after each treatment)*

- This turbine can be immersed (for a maximum of 15 minutes) during the pre-disinfection stage on condition that the following requirements are met:
- Use a detergent-disinfectant specifically designed for pre-disinfection of dynamic dental instruments which contains neither sodium oxide nor aldehyde and which is aluminium compatible.
- Use the concentration recommended by the pre-disinfection product manufacturer.
- Observe the immersion time recommended by the pre-disinfection product manufacturer.
- When the immersion phase is over, rinse thoroughly to drain away the solution/residues from inside the instrument.
- Immediately after rinsing, dry the inside of the instrument meticulously with compressed air.
- Before sterilisation, lubricate the instrument using a "500l" spray can.
- The final three stages of the process (rinsing, drying and lubrication) may be carried out using automated devices available on the market by following the manufacturer's instructions.
- MICRO-MEGA® recommends the use of X-CIDOL® 2 as a pre-disinfection.

#### Automated cleaning

##### *(Steps to take after each treatment)*

- In order to remove the need for handling, use X-Cid® after each patient.
- These devices clean, pre-disinfect and lubricate your instruments without the need for disassembly.
- These devices guarantee fast-acting disinfectant, fungicidal, tuberculocidal and virucidal (HIV, HBV...) effects.
- Regular use of X-Cid® prolongs the useful life of your instruments.

#### Disinfection

Do not use products which contain acetone, chlorine, or bleach as disinfectants! When using disinfectant wipes or jets, follow the manufacturer's instructions.

### 3 Reprocessing Methods

#### Maintenance (Lubrication)

- Insert the "500l" spray can into the back of the turbine.
- Spray repeatedly for a few seconds each time until the liquid expelled is clean.
- Turn the turbine before draining off the excess oil.
- For extended periods of use, repeat these steps every 10 minutes.
- Do not use vinegar to clean the spray nozzles.

#### Packaging

- Once it has been prepared (manually or automatically) for sterilisation and lubrication, place the turbine in an airtight sterilisation bag:
- Individually: A standard packaging material may be used. Ensure that the bag is large enough to hold the instrument without it pulling on the seal.
  - In batches: The instruments may be placed in a special purpose instrument carrier or in separate, general purpose sterilisation trays. Wrap the trays in the appropriate way.

#### Sterilisation

MICRO-MEGA<sup>®</sup> dynamic instruments can be sterilised in a moist heat autoclave for 18 minutes at between 135°C (275°F) and 138°C (280.4°F) with pressure of 2.2 bars inside each bag.

#### Sterilization in a sterilizer compliant with EN 13060

#### Contact corrosion from moisture can occur

Immediately remove the product from the steam sterilizer after the sterilization cycle is complete!

The turbine has a max. temperature resistance of 138°C (280.4°F).

### 3 Reprocessing Methods

#### **Attention:**

Lubricate the turbines before sterilization!

Be aware of the possible risks of the potentially moist environment of the sterilizer when leaving the instruments for extended periods in the chamber with an unopened door. Follow the manufacturer's operating manual.

#### **Storage**

- Observe the manufacturer's recommended shelf life for the sterilisation bags.
- Remove the turbine from the sterilisation bag immediately prior to use.
- Prepared turbines should be stored protected germ-free from dust in a dry, dark and cool room.

#### **Additional information**

When sterilising several instruments in one autoclave cycle, ensure that the autoclave's maximum load is not exceeded (see the instructions for the autoclave used).

- The sterilization bag must be large enough for the turbine so that the bag is not stretched.  
The quality and use of the packaging of the items to be sterilized must satisfy the applicable standards and be appropriate for sterilizing. Shrink-wrap the turbine separately in a sterilizable package!

## 4 Installation and Commissioning

### ! Checking the air pressure:

For operating turbine an operating pressure of minimum 2.1 bar (30 psi) is required. 2.8 bar (40 psi) is recommended. A higher operating pressure is automatically reduced inside the turbine.

The air consumption is 39 - 55 NI/min.

Pressure readings for the turbine:

Drive air = 2.1 - 4.2 bar (30 - 60 psi)

Recommendation = 2.8 bar (40 psi)

Return air < 0.5 bar (7 psi)

Water = 0.8 - 2.5 bar (11 - 36 psi)

Spray air = 1.0 - 4.0 bar (14 - 57 psi)

The minimum amount of water required for spray cooling is 50 ml per minute. An insufficient amount of water may cause overheating of the tooth.

Checking the turbine coupling:

All 5 O-rings must be in place and should not be damaged. Missing or damaged O-rings should be replaced or they may result in premature wear and failure.

### Possible hazard caused by non-sterile products

Danger of infection to the practitioner and patient.

The turbine must be sterilized before first use and after each use.

#### 4 Installation and Commissioning

**Damage caused by contaminated and moist cooling air**

Contaminated and moist cooling air can cause malfunctions and premature bearing wear.  
In general, ensure that a dry, clean and uncontaminated supply of cooling air is provided in accordance with EN ISO 7494-2.

**Burning hazard caused by hot turbine**

When the turbine overheats, burning in the oral area can occur.  
Never contact soft tissue with the turbine!  
The user must use the turbine by knowing that an overheating is possible.

#### 5 Operation

**Connection to units:**

A supply of dry, clean and contamination-free compressed air according to EN ISO 7494-2 must be provided. Contaminated and humid compressed air will result in premature bearing wear. Screw the turbine coupling onto the tubing and tighten firmly.

Changing the turbine:

To change the turbine, place the turbine straight onto the Connex coupling in an axial direction and push until it clicks audibly into place. Also, compatible with KaVo Multiflex® couplers.

Check for secure seating on the coupling by pulling on the turbine. Inexact coupling, especially during the shutoff delay period, may destroy the coupler's high-pressure lamp or reduce its service life.

## 5 Operation

### **Disconnecting:**

Hold the coupling and pull off the turbine while turning it slightly in an axial direction.

### **Diamond and Carbide Use:**

Only use burs that have shafts compliant with ISO 1797-1 type 3, made from steel or tungsten carbide and that meet the following criteria:

- Shaft diameter from 1.59 - 1.60 mm
- Overall length up to max. 25 mm
- Shaft clamping length of at least 11 mm
- Cutting edge diameter up to max. 2 mm
- Short shank burs are recommended for the turbine STATIS SL model 101.1

Ensure that the product is used for its intended purpose only and in compliance with the manufacturer's instructions. Burs that do not meet the above criteria must not be used. Never use burs with worn shafts as these may drop out during treatment, causing injury.

### **Inserting the bur:**

Press the push-button firmly with the thumb. At the same time, insert the bur as far as it will go. Check the function of the chuck by pulling on the bur. Use gloves or another type of finger protection to test, insert and remove in order to avoid any risk of injury and infection.

## 5 Operation

### Removing the bur:

After bur stop, press the push-button firmly and remove the bur. Do not press the push-button while the bur is in motion. Actuation during operation may result in damage to the chuck system and cause injury. It is essential to avoid contact between soft tissue and the head back-cap of the turbine. In the event of turbine damage, heating-up of the head and associated burns are possible.

Remove bur from the turbine when finished with treatment. Placing the instrument in its holder with a bur remaining in the chuck may result in injury and infection. Secure the turbine to prevent it from dropping since this may damage the turbine.

### Damage of the turbine

A turbine that is not connected correctly can become loose from the supply hose and become disconnected. Before each treatment, check that the turbine is safely connected to the supply hose by gentle pulling on the turbine.

### Loss of functionality

Pressing the footswitch while connecting or disconnecting the turbine can cause damage to the turbine and to the supply hose.

Do not connect or disconnect the turbine while the footswitch is operating.

### Spray cooling

#### Risk due to insufficient quantities of water for spray cooling

If the supply of cooling water is insufficient then this can cause the tooth to heat up.

A water supply of at least 50 ml/min is required for the spray cooling.

## 6 Maintenance and Servicing

### **Maintenance**

Only use maintenance products and/or components approved by MICRO-MEGA®.

The use of other products and components can void the warranty.

The turbine is supplied in non-sterile form. Before use, please follow the instructions under this heading.

Check that the sterilizer and the water that is used are clean. After each sterilization cycle, remove the device from the sterilizing apparatus immediately, in order to reduce the risk of corrosion.

**MICRO-MEGA® recommends annual service for all instruments.**

### **Servicing**

Improper servicing and maintenance may result in premature wear and breakdown of the product.

The manufacturer will assume responsibility for proper functioning of these products only if the maintenance products previously named in these instructions are used.

Other cleaning/lubricating products may damage the instruments and void the warranty.

## 7 Specifications

Operating pressure:	2.1 - 4.2 bar	(30 - 60 psi)
Recommendation:	> 2.8 bar	(> 40 psi)
Air consumption:	39 - 55 Nl/min	
Idling speed STAT/S ML:	330,000 - 410,000 rpm	
Idling speed STAT/S SL:	380,000 - 460,000 rpm	
Recommended working pressure:	2 - 3 N	
Weight:	77 g	(0.169 lb)
Dimensions:	120 x 40 x 16 mm	(4.72 x 1.56 x 0.63 inch)
Sound level:	< 70 dBA	

### Ambient conditions

Location:	Permitted in interior rooms
Ambient temperature:	10 to 40°C (50 to 104°F)
Relative humidity:	30 to 75%
Max. altitude:	2,000 m

### Storage and transport conditions

Hazard from starting the turbine after very cold storage. The turbine cannot function.

**Products that are very cold must be warmed to 20°C to 25°C (68°F to 77°F) before use.**

Ambient temperature:	-30 to 70°C	(-22 to 158°F)
Relative humidity:	5 to 95%	
Air pressure:	700 to 1,060 hPa	

Keep dry!

We reserve the right to make technical modifications.

## 8 Guarantee

Your MICRO-MEGA® instrument is an accurate medical device, produced under maximum quality and test conditions. To guarantee personal follow-up and traceability of its medical devices, MICRO-MEGA® would like you to return the guarantee card, duly filled in, within 10 days of your purchase. This return ensures that the procedure for following-up your instrument will be instigated so that MICRO-MEGA® can give you full satisfaction with the service you expect and manage your guarantee correctly.

***You can also register the guarantee for your instrument on your website. [www.micro-mega.com](http://www.micro-mega.com)***

MICRO-MEGA® guarantees your instrument for 24 months against any hidden manufacturing defect, excluding faults due to:

- Improper use.
- Inappropriate maintenance or application, not in accordance with our recommendations and instructions for use.
- Accidental damage (fall, impact, ...).
- Attempted intervention or repair by personnel not approved by MICRO-MEGA®.
- Attempted modification.

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Copyright 2010 MICRO-MEGA®. Tous droits réservés.