# MicroMega K-Files / MicroMega K-Reamers / MicroMega H-Files

**Instructions for use** **EN**

**1. Indication**

The medical devices (K-Files, Hedstrom Files, K-Reamers) are used for catheterisation (01005 to 0110) and root canal shaping (all other diameters) during a non-surgical endodontic treatment. K-Files 008 and 010 are also used in the root canal, initial preparation and root canal permeabilisation.

For use by dental professionals only.

**2. Contraindications**

Apart from children under 2 years of age (ethylene oxide used in the sterilisation process), there are no contraindications to the use for endodontically treating a tooth by orthograde route.

**3. Complications**

In cases of complex canal anatomy, peri-operative risks (instrument breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious process.

**4. Characteristics and warnings**

Ø	Maximum recommended number of uses (if the file is not visually damaged)
< 020	2
≥ 020	5

- Check the condition of the instrument's blade and its fit with the handle before use. If the instrument is damaged or shows signs of wear, do not use.
- Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.

Medical device class according to Directive 93/42/EEC and MDR 2017/745: B.

**5. Clinical claims**

There are no specific clinical claims for the devices other than the sterilisation process; there are no contraindications to the use that corresponds to paragraph 1.

**6. Reprocessing instructions**

- For all metal devices, the use of antiseptic disinfecting and cleaning agents is recommended.
- For your own safety, please wear personal protective equipment (gloves, glasses and mask).
- Do not use cleaning or disinfecting agents using phenol, aldehyde and alkaline compounds.
- Always respect the instructions for use provided by the manufacturer of the products.

**7. Storage and transport conditions**

Place the products in a clean, dry and dust-free environment at a temperature specified by the manufacturer.

**8. Disposal**

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouris) as per good dentistry practices.

**9. Symbols**

STERILEEKO Sterilized using ethylene oxide.

**10. Initial treatment**

Place used products in a container or a wipe with tap water at 20-40 °C and 1.0 % needisher® Septo Active for 5 to 15 min.

Rinse the products with tap water at 20-40 °C for 1 min.

If there is a waiting time before the next step, make sure the device stays moist by placing it in a wet wipe. Do not exceed 1 hour waiting time.

**11. Preparation before cleaning**

If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20-40 °C for at least 1 min with a soft brush until all impurities have been removed.

**12. Visual inspection**

Inspect the used products and discard damaged products (broken, untreated or abnormally bent).

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