
# MicroMega One Glider

Reference: 30006062-A

	Ø	%	L	Stop	Ring	RPM	Maximum Recommended Torque (N.cm)
MicroMega One Glider	16	03	21 - 25 - 31	Yellow	White	250-400 rpm	1.2

## Instruction for use

EN

### 1. Indication

Glide path of root canals during non-surgical endodontic treatment.  
For use by a dental professional only.

### 2. Contraindications

Apart from children under 2 years of age (ethylene oxide used in the sterilization process), there are no contraindications to the use of MicroMega One Glider for endodontically treating a tooth by orthograde route.

### 3. Complications

In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.

### 4. Warnings/precautions

The decision to use an endodontics instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.

• Patients identified as having a risk of infectious endocarditis.

• Contains Nickel and Titanium and should not be used on patients with known allergic sensitivity to these metals.

• Respect the good dental practice in particular the use of a dental dam and gloves.

• Use in continuous rotation at the recommended speed.

• In cases of complex anatomy, the maximum torque may be different from the value recommended by the manufacturer.

• Use according to the recommended procedure (§7).

• Instruments supplied sterile: Check the integrity of the packaging before use. If damaged, do not use the instrument.

• If applicable, respect the expiry date.

• Keep the UDI information which is only on the labelling until the last use.

• Year of manufacture: see labelling.

• Do not use if doubt on the respect of storage conditions.

• Single use instrument: Do not reuse the instrument. If you do, there is a risk of decreased performance and safety characteristics (fracture of the instrument in the tooth).

• In case of doubt concerning the product identification, do not use.

• Check the condition of the instrument before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.

• Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.

### 5. Clinical claims

MicroMega One Glider in normal conditions.

#### Clinical performance:

- Shaping ability: adapted debris removal and cutting efficiency
- Respect of canal anatomy
- Geometry adapted to irrigation
- Maintains the apical foramen in its original position

#### Usability:

- Secure progression and limit the interference

#### Safety:

- Limited risk of breakage
- Limited risk of shaping aberration
- Limited risk of screwing

#### 6. Characteristics

Type 4 rotating instruments under EN ISO 3630-1: 2019 standard:  
• Maximum recommended torque (not required): see table  
'Recommended Torque'.

Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797: 2017 (Type 1).

Number of uses: single use. 4 canals maximum recommended on the same patient

Medical device class according to Directive 93/42/EEC and MDR 2017/745: Ila.

#### 7. Protocol

##### General information on the protocol for use:

After opening the access cavity, sodium hypochlorite must be constantly used and renewed before and after each instrument. The sodium hypochlorite must be placed using a syringe in order to be injected as close as possible to the apical region.

##### Operating procedure :

MicroMega One Glider is introduced into the canal with a slow and regular movement without any apical pressure until the working length is reached.

#### 8. Reprocessing instructions

Not applicable.

#### 9. Storage and transport conditions

#### 10. Disposal

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

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