
MicroMega RECI Glider

Reference: 30006063-A

	Ø	%	L	Stop	Ring	RPM	Angles	Maximum Recommended Torque (N.cm)
MicroMega RECI Glider	16	03	21 - 25 - 31	Red	White	250-400 rpm	CCM 170° CW 60°	4.0

Instruction for use

EN

1. Indication

Glide path of root canals during non-surgical endodontic treatment. Use for a dental professional only.

2. Contraindications

Apart from children under 2 years of age (ethylene oxide used in the sterilization process), there are no contraindications to the use of MicroMega RECI Glider for endodontically treating a tooth by orthograde route.

3. Complications

In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.

4. Warnings/precautions

- The decision to use an endodontics instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
- Patients identified as having a risk of infectious endocarditis.
- Contains Nickel and Titanium and should not be used on patients with known allergic sensitivity to these metals.
- Respect the good dental practice in particular the use of a dental dam and gloves.
- Use in reciprocating rotation at the recommended angles of rotation.
- In cases of complex anatomy, the maximum torque may be different from the value recommended by the manufacturer.
- Use according to the recommended procedure (§7).
- Instruments supplied sterile: Check the integrity of the packaging before use. If damaged, do not use the instrument.
- If applicable, respect the expiry date.
- Keep the UDI information which is only on the labelling until the last use.
- Year of manufacture: see labelling.
- Do not use if doubt on the respect of storage conditions.
- Single use instrument: Do not reuse the instrument. If you do, there is a risk of decreased performance and security characteristics (fracture of the instrument in the tooth).
- In case of doubt concerning the product identification, do not use.
- Check the condition of the instrument before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
- Inform the manufacturer and the national regulatory authority about any serious incident relating to the instrument.

Packaging symbol	STERILE EO		Sterilized using Ethylene Oxide. Do not re-use
------------------	------------	--	------------------------------------------------

5. Clinical claims

MicroMega RECI Glider in normal conditions.

Clinical performance:

- Shaping ability: adapted debris removal and cutting efficiency
- Respect of canal anatomy
- Geometry adapted to irrigation
- Maintains the apical foramen in its original position

Usability:

- Secure progression and limit the interference

Safety:

- Limited risk of breakage
- Limited risk of shaping aberration
- Limited risk of screwing

6. Characteristics

Instrument in reciprocating motion in type 4 under EN ISO 3630-1: 2019 standard: rotation angles are explained in the column <CCW> (Counter-clockwise side) and column <CW> (Clockwise side).

Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797:2017 (Type 1).

Number of uses: single use. 4 canals maximum recommended on the same patient

Medical device class according to Directive 93/42/EEC and MDR 2017/745: IIa.

7. Protocol

General information on the protocol for use:

After opening the access cavity, sodium hypochlorite must be constantly used and renewed before and after each instrument. The sodium hypochlorite must be placed using a syringe in order to be injected as close as possible to the apical region.

Operating procedure :

MicroMega RECI Glider is introduced into the canal with a slow and regular movement without any apical pressure until the working length is reached.

8. Reprocessing instructions

Not applicable.

9. Storage and transport conditions

10. Disposal

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

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