




ROOT CANAL SHAPING

MicroMega

One Curve mini

Rational preparation approach

micromega.coltene.com

 **COLTENE**

- Preserving pericervical dentin
- Removing minimum of tooth structure¹
- Respecting the root canal anatomy²
- Creating optimal space for a deep disinfection and a proper obturation

Protocol MicroMega One Curve mini*



Complete your treatment procedure with matching Guttapercha and Paper Points.

Features

- Small size /taper
- Ø 1 mm wire
- C.Wire heat treatment
- Cutting Efficiency
- Single file, single use
- Sterile blister
- 3 lengths: 21, 25, 31 mm
- Sizes 25/.04, 25/.06, 35/.04, 45/.04



For more information about MicroMega One Curve mini, please scan the QR code.

Order Information

MicroMega One Curve mini, 5 pcs/pack (B5)

60030021	B5 One Curve mini N25 4% L25 Sterile	
60030022	B5 One Curve mini N25 6% L25 Sterile	
60030023	B5 One Curve mini N35 4% L25 Sterile	
60030024	B5 One Curve mini N45 4% L25 Sterile	
60030034	B4 One Curve mini L25 Assorted Sterile	

(1) Survey analysis "Minimally Invasive", Caroline Dort, 2022-01-13

(2) Evaluation of Canal Transportation and Centering Ability of 4 Heat-Treated Nickel-Titanium Systems, JOE, 2020 May;46(5):675-681.doi: 10.1016/j.joen.2020.01.020. Epub 2020 Mar 11., Christina Razcha, Athanasios Zacharopoulos, Dimitris Anestis, Georgios Mikrogeorgis, Giannis Zacharakis, Kleoniki Lyroudia

* Extract from the instructions for use.

© COLTENE – micromega.coltene.com

Coltène/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF/GB
T +44 1444 235486
info.uk@coltene.com

Coltène/Whaledent AG
Feldwiesenstrasse 20
CH-9450 Altstätten
T +41 71 757 5300
info.ch@coltene.com

Micro-Mega SA
12 rue du Tunnel
FR-25000 Besançon
T +33 3 81 54 42 42
commercial@micro-mega.com

